



Qualitative Analysis of Group I Cations

(2015/09/28 revised)


Collect:

- Centrifuge tubes (5)
- Test tube holder
- Dropper
- Latex gloves
- Labels

Prepare:

- Centrifuge (underneath lab bench)
- Test tube rack, test tubes
- * **Concentrated $\text{NH}_3(\text{aq})$: in hood**
- * **Heating the water bath in hood**

Test tube



centrifuge tube





Objective and Techniques


- To learn the techniques of separating and identifying some common cations
- To understand the principles of precipitation and equilibrium of complex formation

Techniques

- Vortex mixer
- Precipitation
- Centrifuge
- Decantation
- Litmus and universal indicator paper



Vortex Mixer



Decantation



Introduction: Qualitative Analysis of Group I~V Cations

Cationic Solutions

- (I) Insoluble chlorides: Hg_2^{2+} , Ag^+ , Pb^{2+}
- (II) Insoluble sulfides in acidic medium: Hg^{2+} , Pb^{2+} , Cu^{2+} , Bi^{3+} , Cd^{2+} , As^{3+} , Sb^{3+} , Sn^{4+}
- (III) Insoluble sulfide or hydroxides in alkaline medium: Al^{3+} , Fe^{3+} , Co^{2+} , Ni^{2+} , Cr^{3+} , Zn^{2+} , Mn^{2+}
- (IV) Insoluble Carbonates: Mg^{2+} , Ca^{2+} , Sr^{2+} , Ba^{2+}
- (V) Soluble cations: NH_4^+ , Na^+ , K^+

Cationic Solution

6 M HCl

Ppt of chlorides of group I cations
AgCl, Hg₂Cl₂, PbCl₂

Solution of group II~V cations

H₂S
pH 0.5

Ppt of sulfides of group 2 cations
HgS, CuS, Bi₂S₃, CdS, PbS,
As₂S₃, Sb₂S₃, SnS₂

Solution of group III~V cations

H₂S
pH 8

Ppt of sulfides or hydroxides of group III cations
Al(OH)₃, Fe(OH)₃, Cr(OH)₃, ZnS, NiS, CoS, MnS

Solution of group IV and V cations

(NH₄)₂CO₃
NH₃/NH₄Cl

Ppt of carbonates of group IV cations
MgCO₃, BaCO₃, SrCO₃, CaCO₃

Solution of group V cations
NH₄⁺, Na⁺, K⁺



Procedure 1-1: Prepare Testing Solutions

Take a centrifuge tube and labeled



Add Hg_2^{2+} , Ag^+ , Pb^{2+}
(2, 2, 3 drops)
to prepare testing soln



* Methods for mixing solutions completely:
(1) Use stirring rod, (2) shake/flick, (3) Use Vortex mixer



Procedure 1-1

Precipitation of Insoluble Chlorides

Cationic solution



step1-1

(1) Add 2 drops 6 M $\text{HCl}_{(aq)}$, stir for 1~2 mins,
(2) Centrifuge and separate supernatant from ppt.

Note: Use Vortex mixer to help mixing thoroughly

Do not add too much $\text{HCl}_{(aq)}$, it may cause the chlorides to dissolve

Ppt 1-1

$\text{AgCl}_{(s)}$, $\text{Hg}_2\text{Cl}_{2(s)}$, $\text{PbCl}_{2(s)}$
(white) (white) (white)

(1) Add 1 d. 6 M HCl with 10 d. of water to wash ppt,
(2) Centrifuge and separate ppt and supernatant

Ppt 1-1

Soln
(Discard)

Soln 1-1

Add 1 drop 6 M $\text{HCl}_{(aq)}$
Check whether precipitation is complete

No

Yes

Repeat steps 1-1

Discard in waste bin



Procedure 1-2 & 1-3

Separation and Identification of Pb^{2+}

Ppt 1-1



Step 1-2

- (1) Add 5 d. distilled water and mix well
- (2) Heat in boiling water bath for several min. in hood
(Water bath should bring to boil to extract PbCl_2 efficiently)
- (3) Centrifuge to separate the ppt. and supernatant

Ppt 1-1 should be extracted with hot water 2~3 times and centrifuged until the extracted solution when added K_2CrO_4 , solution is only slightly yellow

Ppt 1-2

$\text{AgCl}_{(s)}$, $\text{Hg}_2\text{Cl}_{2(s)}$

Step 1-3-2

Extract ppt 1-2 2~3 times and centrifuge to obtain the soln

Ppt 1-2

$\text{AgCl}_{(s)}$, $\text{Hg}_2\text{Cl}_{2(s)}$

Soln 1-2-b

Contains $\text{Pb}^{2+}_{(aq)}$

Soln 1-2-a



Contains $\text{Pb}^{2+}_{(aq)}$

Step 1-3-1

1 d. 6 M HOAc
1 d. 0.5 M K_2CrO_4

Step 1-3-1

**$\text{PbCrO}_{4(s)}$
(Yellow ppt)**






Procedure 1-4

Separation and Identification of Hg_2^{2+}

Ppt 1-2



Step 1-4


(1) Add 2~4 d. 15 M NH_3 (in the hood)

(2) Stir and mix well




(3) Centrifuge to separate ppt. and supernatant



Ppt 1-3




Soln 1-3








Procedure 1-5: Separation and Identification of Ag^+



Add 6 M HNO_3
until solution is acidic



Use stirring glass rod
to drip on litmus paper





Notice

- Prepare hot water bath on hot plate in the hood
- Wear latex gloves throughout experiment
- Take the amount of chemicals according to lab manual, to avoid pollution
- Use test tube rack or test tube holder to transport test tubes or centrifuge tubes; do not use your bare hands
- The experiment produces toxic fumes, **all the heating processes must be done under the hood**
- After centrifuging, solid precipitates and the supernatant should be separated by decantation
- **The liquid waste contains heavy metals and should be collected and discarded into the recycling bin**



Manipulate the Centrifuge

- Use centrifuge tubes in centrifuge, **do not use test tubes**
- Centrifuge tubes should be placed in **opposite sides to keep balanced**
- The lid should be closed during use; the centrifuge should be started from **slow** to check if there are unusual sounds, then the speed can be increased
- If there are unusual sounds or movement in the centrifuge, the power should be turned off first to fix
- There must be at least one person watching the centrifuge when in use
- Centrifuge for approx. **1~2 min**, open lid when the centrifuge has completely stopped