



C26 Recrystallization and Melting Point


(2017/04/24)

Collect:

- Melt-Temp device
- Capillary tube
- Filter paper

Prepare:

- Thermometer
- Capillary tube stand
- Hollow glass tube
- Hot plate
- Erlenmeyer flask (50 mL, 2)
- Hirsch funnel
- Filtering tube
- Rubber stopper
- Water aspirator





C26 Recrystallization and Melting Point

■ Objective:

To learn the techniques of recrystallization and to determine the melting point

■ Technique:

- Electric balance, magnetic stirrer and hot plate
- Recrystallization
- Suction filtration
- Packing capillary tube and m.p. determination

■ Flowchart:

Part I: m.p. determination of pure substances

Part II: recrystallization of benzoic acid

Part II: m.p. determination of crude and purified compd.



Principle

■ Melting point

- The temperature at which it changes from solid to liquid at atmospheric pressure
- When a pure crystalline substance melts, the melting point range should not exceed 1 °C. An impure substance shows a larger melting point range than a pure substance
- The melting point can be used
 - to determine the identity of a known compound
 - to define an unknown compound at later times
 - to determine the purity of a substance



Principle

■ Recrystallization


- Recrystallization consists of dissolving the solid substance in a minimal amount of a suitable hot solvent. Cooling the solution causes the formation of a supersaturated solution, and the pure solid will be crystallized out.
- Ideally, the concentration of the impurity will not reach the saturation level; hence, it will not be crystallized out.
- During recrystallization, the sample adheres to the lattice in an orderly manner, without any impurities included that increases the purity.



Part I Melting Point – Packing Capillary Tube

- Pack 2 capillary tubes for each sample:

(1) benzoic acid (122°C) (2) acetanilide (113°C) (3) benzoic acid/acetanilide (1 : 1)



- Transfer a dry and finely powdered sample onto a piece of weighing paper
- Insert the open-end of a capillary tube into the stack of the sample.



Keep capillary tube open-end up




- Knock the closed-end of the capillary tube on the bench top 3-4 times

- The sample is ca. 2-3 mm height



Part I Melting Point – Determine the Melting Range

- Adjust the voltage to zero, then turn on the power
- Carry out a crude melting point measurement for 3 samples in a faster heating
- * **Start with the lowest m.p. and end with the one with the highest m.p.**
- Cool the Melt-Temp device to lower than the crude m.p. about 15 °C
- Replace the capillary tubes and slowly heat the samples at a rate of 2 °C/min. to determine the m.p.
- Record the m.p. range of the sample, i.e. the temperature range from starting to melt to completely melted



Heating rate

Initial: 10~15 °C/min

Lower than m.p. 10~20 °C: 2 °C/min

Lower than m.p. 2~4°C: 1 °C/min



Part II Recrystallization – Dissolution

- Prepare 2 capillary tubes of crude benzoic acid sample



- Transfer the remaining sample to a 50 mL Erlenmeyer flask after weighing

- Add appropriate amount of hot water to the flask

* **Calculate the amount of hot water needed base on solubility**





- Heat and swirl the solution gently to dissolve
- * **As water is the solvent, heat the solution on hot plate directly instead of water bath**



- Cool soln slowly at room temp for 15 ~ 20 min. to grow crystals



Part II Recrystallization – Crystallization



- Cool the flask in an ice-water bath to decrease the temp. and increase the yield

- Suction filtration to collect crystals
- Wash the crystals with small amount of cold-water
- Suction dry for 10 min.
- Collect the crystals on filter paper and press to dry
- Weigh the dried crystals and determine the yield

- Determine m.p.
(1) crude benzoic acid
(2) recrystallized benzoic acid
(3) **crude acetanilide**



Manipulation of Suction Filtration



- ◆ Fix water-trap bottle and filtering tube with extension clamp
- ◆ Place a rubber stopper between filtering tube and Hirsch funnel
- ◆ Cut the filter paper to cover all of the holes in the bottom of the Hirsch funnel but it must not extend up the sides
- ◆ Moisten the filter paper with small portions of solvent
- ◆ Close the two-way valve, then start and test the suction
- ◆ Start filtering when the filter paper is tightly stuck on the funnel
- ◆ Release the pressure before turn off the water aspirator

- ◆ Fill water in water pump that come in from bottom and flow out from top



Notice

■ Melting point determination

- ◆ The sample inside the capillary tube should be 2~3 mm in height
- ◆ Glass tube for knocking the capillary tube should be kept clean and dry; wash and oven dry after class
- ◆ Record the melting range of sample
- ◆ Replace the capillary tubes in the second run

■ Recrystallization

- ◆ Use the Erlenmeyer flask to dissolve the sample
- ◆ **Use water bath to heat the organic solvents that are flammable**
- ◆ Use the least amount of hot solvent to dissolve the sample to increase the yield
- ◆ Cool the solution slowly at room temp. to grow crystals that increases the purity

- **Recycle:** the benzoic acid, capillary tubes to designated waste bins