



C15 The Electrogravimetric Determination of Copper

Collect

- Platinum gauze and wire
- Stir bar
- Watch glass
- Timer
- DC power supply
- 100 mL beaker



Step 1. Preparation of sample

- Weigh ca. 0.5 g $\text{Cu}(\text{SO}_4)_2 \cdot 5\text{H}_2\text{O}$ into 100 mL beaker
- Add ca. 30 mL of deionized water to dissolve it
- Add another 3 mL of 9 M H_2SO_4 and 2 mL of 6 M HNO_3



Step 2. Preparation of electrodes


- Immerse the Pt electrodes in 6 M HNO₃ for 5 min.
(Operate in fume hood)
- Wash thoroughly with DI water and rinse with ethanol
- Place the electrodes on watch glass and oven-dry at 110 °C for 5 min.
- Cool in desiccators and weigh to the nearest 0.1 mg

NOTE : Pt electrode is very expensive. Please handle with care.





Step 3. Set up the apparatus

(+) Electrode



(-) Electrode



Stir bar

NOTE : Keep the two electrodes away from each other



Step 4. Electrolysis

DC Power supply



- Turn the power on
- Adjust the voltage and current alternatively
- Adjust the current to 1 A
- Turn the magnetic stirring on
- When the blue color of Cu(II) has disappeared, add some distilled water.
- Electrolyze for 15 min at a current of 0.5 A
- Repeat until no further deposition of Cu occurs



Step 5. Determine the mass of Cu

- Wash the Pt gauze electrode with a squirt bottle without turning off the power
- Turn off the power after washing
- Wash the electrode with ethanol and dry with hair-blower
- Cool the electrode and weigh to nearest 0.1 mg
- Calculate the mass% of Cu in the sample
- Immerse the Pt gauze in 6 M HNO_3 to dissolve the Cu and wash with water
- Turn in the Pt electrodes to TA.